

# SALIVA COLLECTION KIT

## Instruction For Use

### INTENDED USE

Used for saliva collection.

### PRINCIPLE

The ergonomically designed collector is easy for saliva collection.

The preservation solution is Hank's Balanced Salt Solution, which can help to stabilize the virus. Phenol red was added to the preservation solution to indicate the preservation state. Antibiotics were added to prevent the growth of bacteria and fungus.

### STORAGE CONDITION AND EXPIRATION

Expiration: view on label

Storage condition: room temperature (2–30 °C).

### KIT CONTENT

1. Saliva Collector
2. Tube pre-filled with 2 ml buffer
3. Barcode stickers
4. Biosafety Bag
5. Test request form
6. Transport bag for sample delivery at the sampling point
7. Instructions for handling of the collected TEST&GO saliva sample
8. Instructions for Use

### USAGE METHOD

1. Rinse your mouth with water before sampling, and do not eat, smoke or brush your teeth before sampling at least 1 hour.
2. Relax your cheeks, spit into the saliva collector, up to 2 ml marker, make sure the amount do not include bubbles.
3. Add the preservation solution pre-filled in tube.
4. Remove the saliva collector and screw the cap.
5. Invert tube for about 10 times to mix.
6. Fill in the necessary information in the test request form and stick the barcode sticker horizontally on the test tube and the marked place in the test request form.
7. Then follow the instructions given in the document "Instructions for handling of the collected TEST&GO saliva sample".



1. Take out the collector and preservation solution.



2. Relax your cheeks, spit into the saliva collector, up to 2 ml marker, make sure the amount do not include bubbles.



3. Add the preservation solution pre-filled in tube.



4. Remove the saliva collector and screw the cap.



5. Invert tube for about 10 times to mix.

### PRECAUTIONS

1. Do not touch the inside of the collector with hands or other items when using the collector.
2. The product is disposal, do not reuse it.
3. Do not use the product when its package is broken.
4. Do not use the product when it expires.
5. The product should be stored in a dry, ventilated, radiation-free, clean room.
6. Do not drink, touch the preservation solution. The solution can be harmful if ingested and may cause irritation if exposed to the skin and eyes.
7. If possible, please wear disposable gloves during whole operation.

**Manufacturing date and expiration date: view on label**

*Breach of the above instructions may affect the result.*



Call centre 0800 800 234  
Mon – Fri 7–17  
[www.synlab.sk](http://www.synlab.sk)

